

ID	001		
Name	Verify functionality for Hardware tab		
Priority	High		
Pre-conditions	Hardware connected, all connections are correct and power is on. Breeze is open.		
Test step	Expected result	Actual Result	Pass/Fail
1. In Breeze, go to <i>Settings</i>	<i>Settings</i> opens and <i>Focus</i> is selected		
2. Select <i>Hardware</i>	<i>Hardware</i> opens		
3. Select Camera in drop-down list <i>Selected Camera</i>	X camera is selected in list		
4. Select camera type in drop-down list <i>Type</i>	X camera type is selected in list		
5. Press <i>Connect</i>	X camera is connected and <i>Information</i> presents information regarding x camera		
6. Select Sample Mover in drop-down list <i>Selected Sample Mover</i>	X sample mover is selected in list		
7. Press <i>Connect</i>	X sample mover is connected		

**Comments**

<b>ID</b>	<b>002</b>		
<b>Name</b>	Verify functionality for Sensor tab		
<b>Priority</b>	Medium		
<b>Pre-conditions</b>	TC001 passed		
<b>Test step</b>	<b>Expected result</b>	<b>Actual Result</b>	<b>Pass/Fail</b>
1. Select <i>Sensor</i>	<i>Sensor</i> is selected		
2. Press <i>Camera</i> button	Sensor image is presenting image with dead pixels		

Comments

<b>ID</b>	<b>003</b>		
<b>Name</b>	Verify functionality for Measurement tab		
<b>Priority</b>	Medium		
<b>Pre-conditions</b>	TC001 passed		
<b>Test step</b>	<b>Expected result</b>	<b>Actual Result</b>	<b>Pass/Fail</b>
1. Select <i>Measurement</i>	<i>Measurement</i> is selected		
2. Set <i>Measurement Length</i> to x	Measurement Length is set		

<b>Comments</b>

<b>ID</b>	<b>004</b>		
<b>Name</b>	Verify functionality for Integration Time tab		
<b>Priority</b>	Medium		
<b>Pre-conditions</b>	TC001 passed		
<b>Test step</b>	<b>Expected result</b>	<b>Actual Result</b>	<b>Pass/Fail</b>
1. Select <i>Integration Time</i>	<i>Integration Time</i> is selected		
2. Press <i>Camera</i> button	Plot and Saturation Ratio is presented		
3. Press <i>White Reference</i> button	White reference is presented		
4. Set <i>Integration Time</i> to X value	<i>Integration Time</i> is set		
5. Check Min and max	Check box is checked and min/max is presented in plot		
6. Use slider below plot to adjust <i>Position</i>	Position is adjusted in plot		

Comments

<b>ID</b>	<b>005</b>		
<b>Name</b>	Verify functionality for Focus tab		
<b>Priority</b>	Medium		
<b>Pre-conditions</b>	TC001 passed		
<b>Test step</b>	<b>Expected result</b>	<b>Actual Result</b>	<b>Pass/Fail</b>
1. Select <i>Focus</i>	<i>Focus</i> is selected		
2. Press <i>Camera</i> button	Image is presented		
3. Set <i>Field of View</i> to x mm	<i>Field of View</i> is set		
4. Press <i>Focus Grid</i> if you have focus grid, if not place focus grid on sample mover and then press <i>Focus Grid</i>	Focus grid image is presented		
5. Check <i>Show line plot</i>	Line plot is presented		
6. Press <i>Scan</i>	Sample Mover is moving and image is presented		
7. Use slider below plot to set <i>Band</i> to x value	<i>Band</i> is set, line plot and focus grid image is adjusted		

**Comments**