

Breeze - Import data guide

Description of options:

- **Use of one dark and white reference (Optional)**

Import one or several image data files and apply only one dark and reference file to all imported files. See *instructions on next page*

- **Camera max signal (Optional)**

The camera maximum signal value. This will be used to for example for example white reference values for saturations.

- **Field of view (Optional)**


The width of the image data files in millimeter. The field of view will be used to converting pixel values into millimeter in for example spatial descriptors

Supported file formats:


File format	Extensions
Breeze xml file	xml
ENVI Files	bil, bip, bsq, raw, img, bin, dat
Matlab files	mat
HDF files	h5, hdf
Image files	jpg, jpeg, wbmp, png, bmp, pbm, pgm, ppm, pcx, tif, tiff, gif, bmp, jp2
SAC file	sac
HIPS File	hips

Import using one dark and white reference

1. **Start Breeze with the shortcut created after installation.**
2. Enter the Record view by **pressing the “Record” button**
3. You will then see the following view. **Select the “Import data” tab.**
Press OK to start Import data wizard



4. **Browse to folder with files to import**
Select one or several files



Press Next


5. Click on **Advanced** to expand the advanced section



Select **“User the same reference for all imported files”**

Browse or enter path for **“Dark reference”**

Browse or enter path for **“White reference”**



Press Finish

6. Your data will be imported into Evince.
If **“Automatically apply workflow to imported files” is selected** then the sample model(s) specified in the **“Analysis Tree”** will be applied to the imported image files.

